

DiscoveryUltra


NEW !!!

**EXCLUSIVELY THROUGH
VIDEO DENTAL CONCEPTS!**

The Best Wireless Intraoral Camera
at the Lowest Price!!

Combine it with the appropriate receiver to capture and save
Images to computers, Digital and Analog Video Displays.

Perfect for **RF SYSTEM^{lab} Einstein** users
and **RF AMERICA DiscoveryCam** users.



Macro Lens



DiscoveryULTRA
Universal USB, RCA,
VGA



DiscoveryULTRA
Universal RCA, VGA,
SD



Channel Selector

ULTRA intraoral camera specifications

Image Sensor	1/4" SONY CCD
Image Resolution	2.0 Megapixels
Output Jack	Varies with Receiver
Lamp	6 Pcs. White LED (5600K)
Depth of Field	3~50mm
Visual Angle	105°
Signal Format	Varies with Receiver
Power Supply	1600mAh Rechargeable Lithium Battery

800-323-2690

VideoDental
CONCEPTS

www.videodental.com