

QuickCam

Direct USB

Advanced Features include:

HIGH END FEATURES AT AN AFFORDABLE PRICE

CAPTURE METHOD

- Dual-sided ambidextrous buttons
- Optional foot-pedal available
- Capture/Save function is programmable in most popular imaging software
- Direct high-speed 2.0 USB connection
- 3m cable with metal connectors

OPTICS

- Optical Train Consisting of 6 elements in 3 groups
- 90 Degree field of view combined with a 4 step dial, renders all popular face smile intra oral and macro modes



ILLUMINATION

- Convenient ON/OFF Feature
- 9 Extremely bright white LED's

QuickCam

Direct USB

PORTABLE HAND-PIECE

- 2.3 OUNCES
- NO DOCKING STATION REQUIRED
- SEAMLESSLY INTERGRATES WITH MOST IMAGING SOFTWARE
- CONVENIENT CABLE/CRADLE KITS AVAILABLE



IMAGE QUALITY

- SONY HIGH RESOLUTION CCD 1/4"
- RESOLUTION 811 X 508

SOFTWARE INTERFACE

- VIDEO FOR WINDOWS, WDM
- OS CAPABILITY
XP/TESTING WINDOWS 7



VideoDental
CONCEPTS

DENTERPRISE INTERNATIONAL INC,
110 EAST GRANADA BLVD. SUITE 207
ORMOND BEACH, FL 32176
TOLL FREE: 800-323-2690
FAX: 386-672-4402
EMAIL: SALES@VIDEODENTAL.COM
WWW.VIDEODENTAL.COM