

SHEATHS AND CAMERA COMPATIBILITY

CAMERA	COMPANY	Population	Sheath 1'	VDC #	Alt Sheath	VDC #	Mfgr #'	Note #
AcuCam (variegated FO tube)	New Image	Original	B50	50-100			20816	
Accucam(old)&Polo ('94-97)	New Image	5000	B50	50-100			20816	
Accucam Concepts 3 ('98-01)	Gendex	2000	B123	50-105			20893	
Accucam Concept IV	Gendex	10000	B116	50-200	B155		20840	
AcuCam Plus	New Image		B50	50-100			20816	
Advantage II	IPS		HardPlastic	Prop				4
American Dental Camera	Cygnus COS-200		B81	50-125			20821	
Aurora	Digital Doc		B184	50-226			20911	
Cad/Cam	IPS/Med Dynamics		HardPlastic	Prop				4
CamaLot (aka OraCam)	Flow Xray/Sometech		B193	50-227	B173	50-130	20915	
Camera Caps (Xray viewing)	Flow Xray		#16810	50-240				
Camrex 191 (aka MiniCam+)	DentAmerica		B173	50-130			20909	
Camrex 191A, 192	DentAmerica		B195	50-229			20916	
Celebrity	Healthco/WA	1000	B40	50-115			20812	1, 2
ClearCam (see QuickCam)								
CompuCam (see QuickCam)								
Copernicus	RF System Labs		B180	50-136			20964	
Cudaprism	Cuda	100	B73					
Cygnascope	Zila/IDT	1000	B81	50-125			20821	
DP-6/DP-15	RF System Labs		B117	50-120				
DP-6 ver.2	RF System Labs		RF Short	50-137	B180	50-136	20964	
DentaCam	Fuji	1000	B-179	50-135			20963	
DexCam (see EasyDoc)								
DiamondCam (see EasyDoc)								
Digital Doc	Digital Doc		B151	50-225			20883	
Dr. Camscope	Mikky Intl	100	B168					
DX Small #1 Sensor	Various	1000	B87	50-100			20824	
DX Large #2 Sensor	Various	1000	B88	50-113	B79		20825	
EasyCam	Schein/ICS	5000	B40	50-115			20812	2
Easycam Ultra Lite II	ICS		B148	50-180			20884	2
EasyDoc	ICS		B117	50-120			20831	
EasyDoc II	ICS		B148	50-180			20884	
EasyOp (see EasyDoc)								
EasyVu (see EasyDoc)								
Einstein (Satellite Scope)	RF System Labs		RF Short	50-137	B180	50-136	20964	
Einstein - Spark	RF System Labs		RF Short	50-137	B180	50-136	20964	
EuroDoc (see EasyDoc)								
Explorer	Expanded Optics	500	B73					
F-5W	RF System Labs		RF Short	50-137	B180	50-	20964	
Insight	Dentsply/New Image	7000	B83	50-145			20856	
Intracam	ICS	5000	B40	50-115			20812	2
Litex 695 (Camrex 191A)	DentAmerica		B173	50-130			20909	
LogiCam (see QuickCam)								
MedRx (see QuickCam)								
Miharu (v1) Home Use Cam))	RF System Labs		B86	50-150			20855	
Miharu (v2) Production Vers))	RF System Labs		RF Short	50-137	B180	50-136	20964	
MiniCam / MiniCam USB	VDC	1500	B173	50-130			20909	
MiniCam Plus (Camrex)	VDC / DentAmerica		B173	50-130			20909	
MiniCam Ultra(aka Sopro 595)	VDC	350	B86	50-150			20855	
NetCam+	IPS		HardPlastic	Prop				4
Omni	Air Techniques	1000	B84		B73			
OraCam	VDC / Sometch	500	B193	50-227	B173	50-130	20915	
OralScan	LesterDine/VDS	1000	0°		90°		Prop.	4, 5
OralVideoScope (OVS)	VDC	1000	0°		90°		Prop.	4, 5
OralVision	Zila/IDT	2000	B53	50-170	B164			
OralVision 1000	Cygnus		B179	50-135			20963	
Perspective	Dentsply/NewImage	1000	B40	50-115			20812	2

Prism	Prodentec/ICS	500	B148	50-180			20884	2
ProCam	Prodentec/ICS	500	B148	50-180			20884	2
ProScope 1000	DentX/Prodent	1000	B117	50-120			20831	
Provision	DIS	500	B53					
QuickCam (DDS/555/575)	VDC / Sopro	200	B116	50-200			20840	
QuickCam PC / "2001"	VDC / Owandy	100	B179	50-135			20963	
QuickRay Sensor	VDC		B87	50-110			20824	1
Reveal	Patterson/WA	2000	B40	50-115			20812	2
SS-21 C-Line Satellite Scope	RF System Labs		RF Short	50-137				
SS-24 Satellite Scope	RF System Labs		RF Short	50-137				
Schick 1	Schick Technologies	500	B160					
Schick 2	Schick Technologies		B126		B123	50-105	20893	
Smilie575 (see QuickCam)								
Spectravue	ICT	3000	Metal	Prop.	B100	50-215		4
Sopro 555 (see QuickCam)								
Sopro 575 (see QuickCam)								
Sopro 595 (see MiniCamUltra)								
SuniCam USB	Suni Imaging		B201				20969	
Suni Sensor #0	Suni Imaging		B87	50-110			20824	
Suni Sensor #1	Suni Imaging		B87	50-110			20824	
Suni Sensor #2	Suni Imaging		B131				20861	
SuperCam (see EasyDoc)								
TeliCam	DMD	5000	B108	50-220			20810	
TeliCam Elite (Blu/Wht Case)	DMD	5000	B108	50-220			20810	
Toothpix	VDC	500	B116	50-200			20840	
Trophy	Trophy	1000	B58		B96			
TrueVision (see Advantage II)								
UltraCam / UltraCam II	ADT/Ultracam	2000	HardPlastic	Prop		50-230		4, 6
VidaCam	ICS	5000	B40	50-115			20812	2
VideoScope2000 (V2000)	VDC	500	0°		90°		Prop.	4, 5
ViperCam	Integra/Dynamic	2000	B116	50-200			20840	
ViperCam II	Integra/Dynamic		B172	Prop.				4
VistaCam	Air Techniques	3000	B84	B180	B73			
Viola	Dynamic Dental		B116	50-200			20840	

Notes:

1. WA = Welch Allyn
2. ICS Order
- 3
4. Prop. = Proprietary
5. Metal Rotator by Cuda Products
6. Sheath Mfgr Medical Dynamics